

# Vaccine Storage Safety !

Be sure vaccines protect.  
Handle with care!

## Vaccine Refrigerators

Refrigerators used to store vaccines must be selected carefully and used properly.

Any refrigerator used for vaccine storage should be:

✓	Able to maintain recommended vaccine storage temperatures (between +2° C to +8° C) through all seasons
✓	Large enough to hold the year's highest monthly inventory (including influenza season)
✓	Equipped with a calibrated thermometer or data logger inside each storage compartment
✓	Dedicated to the storage of vaccines only
✓	Placed in a secure location away from unauthorized and public access

**Bar-style fridges** (small, single door) are unpredictable for maintaining temperature and should **not be used** for vaccine storage.